

Medela Breast Pumps

The #1 recommended Breast Pump Brand*



Advanced Double Breast Pump



Pump in Style® Advanced Starter Set



Pump in Style® Advanced - Backpack & On-the-go Tote



Pump In Style® Advanced The Metro Bag™



Freestyle® Breast Pump



Sonata® Breast Pump



Symphony® Breast Pump

SKU	101035078	101035077	57062 / 57063	57036	101034712	58200	0240108 / 0240208
Description Marketing or Packaging	Medela's Advanced Personal Double Breast Pump delivers an effective pumping solution which includes researched performance of 2-Phase Expression technology, adjustable speed and vacuum and double pumping capabilities.	Medela's Pump in Style Advanced Starter Set delivers effective pumping in a modular solution. Offering trusted and research based performance with 2-Phase Expression in a streamlined format customizable to mom's changing lifestyle.	Medela's award-winning Pump in Style Advanced breast pump designed for moms who need true portability when pumping on-the-go. Containing Medela's trusted research based performance of 2-Phase Expression technology providing mom everything she needs to get started and stay organized while breast pumping.	The same great Pump in Style Advanced Breast Pump with a premium carry bag and extra accessories. This double electric breast pump, designed for portability, helps mom balance breast pumping with her goals & everyday life. It contains Medela's trusted research based performance of 2-Phase Expression technology in an all-in-one breast pump system.	Medela's Freestyle Breast Pump with its compact and lightweight system along with rechargeable battery and belt clip is designed to make breastfeeding easier for true mobility. Offering trusted performance of 2-Phase Expression technology, Freestyle provides mom the ultimate double pumping experience by helping her maintain her active lifestyle.	Medela's Sonata Breast Pump offers our first connected solution to empower and support moms, breast milk feeding success through smart connectivity for confidence and convenience, innovative performance features designed to ensure optimal efficiency in every session and quiet operation for a comfortable experience.	Medela's Gold-Standard hospital-grade (multi-user) breast pump designed for breast milk initiation, building and maintenance.
Technology							
Daily Use	●	●	●	●	●	●	●
Variable Vacuum/ Speed Settings	●	●	●	●	●	●	●
mmHg - Vacuum	Combined & Fully Adjustable -20 - -270 mmHg	Combined & Fully Adjustable -20 - -270 mmHg	Combined & Fully Adjustable -20 - -270 mmHg	Combined & Fully Adjustable -20 - -270 mmHg	Combined & Fully Adjustable -20 - -270 mmHg	Combined & Fully Adjustable -50 - -250 mmHg	Combined & Fully Adjustable -30 - -270 mmHg
Cycles Per Minute	54-120 cpm	54-120 cpm	54-120 cpm	54-120 cpm	45 - 120 cpm	54-120 cpm	54-120 cpm
Carry Bag	Motor bag only (carry bag sold separately)	Motor bag only (carry bag sold separately)	All-in-one portable black microfiber shoulder bag or backpack.	All-in-one portable black microfiber shoulder bag includes work surface and mesh bags for part organization.	All-in-one portable black microfiber shoulder bag	All-in-one portable black microfiber shoulder bag	Hard protective carry case
Removable Cooler & Ice Pack	Sold Separately	Sold Separately	●	●	●	●	Sold Separately
Breast Milk Containers	2 (5 oz)	2 (5 oz)	4 (5 oz)	4 (5 oz)	4 (5 oz)	4 (5 oz)	2 (5 oz) included in the double pumping kit
Battery	Sold Separately (8AA external battery pack)	8AA external battery pack (AA not included)	8AA external battery pack (AA not included)	8AA external battery pack (AA not included)	Rechargeable Lithium-Ion 3 hrs.	Rechargeable Lithium-Ion 1 hr.	Symphony Breast Pump with Rechargeable Battery 0240208 1 hr.
Digital Display	--	--	--	--	●	●	●
Timer	--	--	--	--	--	Count-up and Count-down	--
Hospital-Grade (multi-user)	--	--	--	--	--	--	●
Hospital Performance	--	--	--	--	--	Clinically Tested Hospital Performance Pumping Rhythms designed for maintaining milk supply	Clinically Tested Hospital Performance Pumping Rhythms designed for initiation, building & maintaining milk supply
Overflow Protection	At the Pump	At the Pump	At the Pump	At the Pump	At the Kit	At the Kit	At the Pump
Memory Feature	--	--	--	--	●	--	--
Bluetooth® Connectivity	--	--	--	--	--	●	--

Content appearance may vary from picture.

*QuintilesIMS ProVoice Survey: Nov-Oct 2017. The Bluetooth word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Medela is under license. Medela wordmark and logo, Pump in Style, Freestyle, Sonata, Symphony, and 2-Phase Expression are registered in the U.S. Patent and Trademark Office. 1549080 B 0418 © 2018 Medela 4551067